Biochemistry Lab, Clinical Biochemistry Lab

Document Name Directory of services

SJH/CBU/Clinical Biochem Lab/DOS

CLINICAL BIOCHEMISTRY LAB, OLD EMERGENCY BUILDING, $3^{\rm RD}$ FLOOR, VMMC AND SAFDARJUNG HOSPITAL, NEW DELHI

Document Number

OPENING HOURS: 8:30 AM TO 4 PM

		DIREC	TORY OF SE	RVICES		
S. No.	Name of the Test	SPECIMEN CRITERIA Amount Type of		Instructions	NABL Accreditation	
1.	Glucose F	2 ml blood	Vacutainer Fluoride (Grey capped)	8-10 hrs of fasting	Yes	
2.	Glucose PP	2 ml blood	Fluoride (Grey capped)	2 hrs after taking breakfast	Yes	
3.	Glucose R	2 ml blood	Fluoride (Grey capped)	Random sample	Yes	
4.	OGTT	2 ml blood	Fluoride (Grey capped)	2 hrs after taking 75 gms glucose	Yes	
5.	Urea	2 ml blood	Plain (Red capped)/ Gel (Yellow capped)	Random sample	Yes	
6.	Creatinine	2 ml blood	Plain (Red capped)/ Gel (Yellow capped)	Random sample	Yes	
7.	Uric Acid	2 ml blood	Plain (Red capped)/ Gel (Yellow capped)	Random sample	Yes	
8.	Calcium	2 ml blood	Plain (Red capped)/ Gel (Yellow capped)	Random sample	Yes	
9.	Phosphorus	2 ml blood	Plain (Red capped)/ Gel (Yellow capped)	Random sample	Yes ONTROLLED DOCU	ME

Issue No.: 01 Issue Date:		Amendment Date:	CONTROLLED DOCUMEN	
07 04 2021		23 01 2023	g g	
Prepared by: Dr.	Reviewed by: Dr Vand	ana Approved by:	Dr Anita Rani	Page no. 1
Khageshwar Mohado	Saini Volume	<u>/</u> /	4	

Biochemistry Lab, Clinical Biochemistry Lab

Document Name
Directory of services

SJH/CBU/Clinical Biochem Lab/DOS

10.	Sodium	2 ml blood	Plain (Red capped)/ Gel (Yellow capped)	Random sample	Yes
11.	Potassium	2 ml blood	Plain (Red capped)/ Gel (Yellow capped)	Random sample	Yes
12.	Chloride	2 ml blood	Plain (Red capped)/ Gel (Yellow capped)	Random sample	Yes
13.	Bilurubin, Total	2 ml blood	Plain (Red capped)/ Gel (Yellow capped)	Random sample	Yes
14.	Bilurubin, Direct	2 ml blood	Plain (Red capped)/ Gel (Yellow capped)	Random sample	No
15.	Total Protein	2 ml blood	Plain (Red capped)/ Gel (Yellow capped)	Random sample	Yes
16.	Albumin	2 ml blood	Plain (Red capped)/ Gel (Yellow capped)	Random sample	Yes
17.	AST	2 ml blood	Plain (Red capped)/ Gel (Yellow capped)	Random sample	Yes
18.	ALT	2 ml blood	Plain (Red capped)/ Gel (Yellow capped)	Random sample	Yes

Issue No.: 01	Issue Date: 01/04/2021	Amendment No: (3)	Amendment Date: 23 01 2023	CONTROLLE	D DOCUMENT
Prepared by: D Khageshwar M		Reviewed by: Dr-Vand Saini	dana Approved by:	Dr Anita Rani	Page no. 2

Biochemistry Lab, Clinical Biochemistry Lab

Document Name
Document Number

Document Number

Document Number

Directory of services

SJH/CBU/Clinical Biochem Lab/DOS

19.	ALP	2 ml blood	Plain (Red capped)/ Gel (Yellow capped)	Random sample	Yes
20.	Amylase	2 ml blood	Plain (Red capped)/ Gel (Yellow capped)	Random sample	Yes
21.	CK	2 ml blood	Plain (Red capped)/ Gel (Yellow capped)	Random sample	No
22.	CK MB	2 ml blood	Plain (Red capped)/ Gel (Yellow capped)	Random sample	No
23.	Body Fluid Glucose (Ascitic/ Pleural/ Pericardial /Synovial)	1-2 ml of body fluid	Sterilized container	Random sample	No
24.	Body Fluid Protein (Ascitic/ Pleural/ Pericardial/ Synovial)	1-2 ml of body fluid	Sterilized container	Random sample	No
25.	LDH	2 mL blood	Plain (Red capped)/ Gel (Yellow capped)	Random sample	No
26.	Magnesium	2 mL blood	Plain (Red capped)/ Gel (Yellow capped)	Random sample	Yes
27.	Total Cholesterol	2 mL blood	Plain (Red capped)/ Gel (Yellow	10-12 hrs fasting	Yes TROLLED DOG

	×		(Yellow			CONTRO	DLLED DOCUM	ENT
	7/04/2021	Amendment N	No: Am	1 4	ent Date:	CONTROLLE	D DOCUMENT	
Prepared by: Dr. Khageshwar Mohado		Reviewed by: Saini	Dr Vandana		proved by:	Dr Anita Rani	Page no. 3	

Biochemistry Lab, Clinical Biochemistry Lab

Document Name

Directory of services

SJH/CBU/Clinical Biochem Lab/DOS

			capped)		
28.	HDL Cholesterol	2 mL blood	Plain (Red capped)/ Gel (Yellow capped)	10-12 hrs fasting	Yes
29.	LDL Cholesterol	2 mL blood	Plain (Red capped)/ Gel (Yellow capped)	10-12 hrs fasting	No
30.	Triglyceride	2 mL blood	Plain (Red capped)/ Gel (Yellow capped)	10-12 hrs fasting	Yes
31.	Spot Urinary creatinine	5 mL of urine	Urine container	Random sample	No
32.	Spot Urinary glucose	5 mL of urine	Urine container	Random sample	No
33.	24 hrs Urinary creatinine	Total amount collected in 24hrs	Urine container	Strictly timed 24hrs collected sample	No
34.	24 hrs Urinary calcium	Total amount collected in 24hrs	Urine container	Strictly timed 24hrs collected sample	No
35.	24 hrs Urinary sodium	Total amount collected in 24hrs	Urine container	Strictly timed 24hrs collected sample	No
36.	24 hrs Urinary potassium	Total amount collected in 24hrs	Urine container	Strictly timed 24hrs collected sample	No
37.	24 hrs Urinary chloride	Total amount collected in 24hrs	Urine container	Strictly timed 24hrs collected sample	No

Issue Date: Amendment No. a. Amendment Date: CONTROLLED DOCUMENT

07/04/2021		endment Date:	CONTROLLE	D DOCUMENT
Prepared by: Dr. Khageshwar Mohado	Reviewed by: Dr Vandana Saini	Approved by:	Dr Anita Rani	Page no. 4

VardhmanMahavir Medical College &Safdarjung Hospital Biochemistry Lab, Clinical Biochemistry Lab Document Name Directory of services SJH/CBU/Clinical Biochem Lab/DOS

38.	24 hrs Urinary phosphate	Total amount collected in 24hrs	Urine container	Strictly timed 24hrs collected sample	No
39.	24 hrs Urinary urea	Total amount collected in 24hrs	Urine container	Strictly timed 24hrs collected sample	No
40.	24 hrs Urinary uric acid	Total amount collected in 24hrs	Urine container	Strictly timed 24hrs collected sample	No
41.	24 hrs Urinary glucose	Total amount collected in 24hrs	Urine container	Strictly timed 24hrs collected sample	No

TURN AROUND TIME, CLINICAL BIOCHEMISTRY LAB, VMMC AND SAFDARJUNG HOSPITAL, NEW DELHI

S. NO.	ZONE	TURN AROUND TIME
1.	OPD	1 Day
2.	Ward/IPD	7 hours
3.	ICU	6 hours
4.	Urgent samples	2.5 hours

CONTROLLED DOCUMENT

Issue No.: 01 Issue	Date: Amend	Iment No: 🌣 📗 Am	The second of th		CONTROLLED DOCUMENT	
07/0	4 2021	2	3 01 2	0 23		
Prepared by: Dr.	Review	ved by: Dr-Vandana	Approve	ed by: Dr Anit	a Rani	Page no. 5
Khageshwar Mohado	Saini	1 Joursey		4	-	
		0, 0		011/		